

The Dream Team...

of foreign body management devices



us endoscopy

listening...and delivering solutions®

The Dream Team...

For over 20 years, US Endoscopy has led the innovation and development of the specialty tools needed for successful Foreign Body retrieval. We are proud to offer a complete and extensive portfolio of retrieval solutions to ensure you have the right product for all clinical situations.



Guardus® overtube

Gentle protection and enhanced safety during foreign body extraction.

BX00711146 (25cm length)
BX00711148 (50cm length)
standard scopes (8.6-10.0mm scope OD)

BX00711147 (25cm length)
BX00711149 (50cm length)
therapeutic scopes (9.8-11.7mm scope OD)



Raptor® grasping device

Offers the versatility needed for retrieving a variety of foreign bodies and stents.

BX00711177
sheath diameter 2.4mm
length 230cm

BX00711178
sheath diameter 1.9mm
length 200cm



Talon® grasping device

Grips into the food impaction for a swift retrieval.

BX00711175
sheath diameter 2.5mm
length 160cm

Roth Net® retrievers

US Endoscopy offers a broad portfolio of 20 Roth Net® retrievers to support a variety of clinical applications and retrieval needs. From foreign bodies to challenging food bolus cases, the following will cover you from top to bottom.

NEW Roth Net® retriever - select

BX00711169* / BX00711170
net size 3.0cm x 6.0cm
sheath diameter 2.5mm
length 230cm

Roth Net® retriever - standard

BX00711050 / BX00711052
BX00711162*
net size 3.0cm x 6.0cm
sheath diameter 2.5mm
length 230cm

Roth Net® retriever - mini

BX00711057 / BX00711168*
net size 2.0cm x 4.5cm
sheath diameter 1.8mm
length 160cm

Roth Net® Platinum® retriever - universal

BX00715050 / BX00715060*
net size 4.0cm x 5.5cm
sheath diameter 2.5mm
length 230cm




Roth Net® Platinum® retriever - food bolus

BX00711155
net size 4.0cm x 5.5cm
sheath diameter 2.5mm
length 160cm

Roth Net® retriever - maxi

BX00711059
net size 4.0cm x 8.0cm
sheath diameter 3.0mm
length 160cm

BX00711165*
net size 4.0cm x 8.0cm
sheath diameter 3.0mm
length 230cm



*sterile

Falcon® rotatable retrieval basket

A great utility tool to have on hand for removal of long objects and other foreign bodies not suitable for the Roth Net® retriever.

BX00711141
basket size 3.5cm x 6.0cm
sheath diameter 2.5mm
length 230cm



US Endoscopy offers Continuing Education and in-servicing
to ensure you and your staff are prepared for any clinical
challenge – contact us today to schedule.

www.usendoscopy.com/education



US Endoscopy is a subsidiary of STERIS Corporation.
©2017 US Endoscopy. All Rights Reserved.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design,
and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226
email info@usendoscopy.com

www.usendoscopy.com

761372 Rev. A